

体外诊断试剂市场调研申请表

拟了解项目	一、适用于贝克曼 Navios 流式细胞仪的流式检测抗体： CD2,CD4,CD5,CD7,CD8,CD10,CD13,CD14,CD16,CD19,CD20,CD33,CD34,CD38,CD59,CD71,CD11b,HLA-DR,MPO,TdT,CD61,CD9,CD15,CD23,CD41,CD43,CD58,CD103,IgM,TCR$\gamma\delta$,Kappa,FMC7,CD3,CD29,CD62L,CD62P,CD45RO,CD49B,CD49E,CD127,CD64,CD45RA,CD1a,CD57,CD303,TCR$\nu\delta$,IFN-γ,IgG1,TRBC1,CD40,CD22,CD24,CD56,CD117,CD235a,CD79a,CD158e, TCR$\alpha\beta$,ZAP-70,CD123,CD11C,CD25,CD55,CD42B,Lambda,CD28,IL-4,TIA-1,BCL-2,BCL-6,Ki67,CD30,CD66c,CD26,CD203c,CD138,CD99,CD304,CD105,CD94,CD69,CD45,CD158a,7AAD, CD27 ,CD158b,CD158i,CD200,CD36,CD159a,CD81,CD279,FLAER,CD38 Nanobody, CD161,CD159c,CD300e,CD180,CD148,NG2,CD66c,GD2,CD326, TCR$\nu d1$, TCR$\nu d2$, CD274, CD95, BCMA,CK,CD141,CD1c,Lin, clec9, CD133, ROR1, Granzyme B,perforin, IgD,CD197,TCRVβ试剂盒 二、其他相关流式试剂 破膜剂, CD34 干细胞计数质控物, 流式细胞溶血剂, 流式细胞鞘液。 三、适用于贝克曼 Navios 流式细胞仪的耗材： 流式上样管, FLOW CHECK PRO, FLOW SET PRO, 流式细胞仪
-------	--

	clenz 清洗液。
该项目的意义	用于血液系统疾病流式细胞免疫分型,辅助白血病、淋巴瘤、多发性骨髓瘤等疾病的诊断与鉴别诊断。